

# SWERECO



**BELMONT  
MEDICAL  
TECHNOLOGIES**

Saving Lives. Together.

## CRITICOOL®

Safe, Precise, and Effective<sup>1,2</sup>

PATIENT TEMPERATURE MANAGEMENT



CRITICOOL™ ANVÄNDS IDAG INOM  
FÖLJANDE TERAPIER/OMRÅDEN

- Cardiac arrest, intensivens
- Asfyxi, neonatalintensiven
- TBI patienter med mild hypotermi, neurointensiven
- Operation
- Sepsis

# SWERECO

## Orderinformation

Art. nr.	Beskrivning	Förp.
Criticool	CritiCool™	1 st
200-00147	Slang till CureWrap™, 2-by-3	1 st
200-00109	Slang till ThermoWrap™, 2-by-2	1 st
014-00028	Kärnadapter (för engångssensorer)	1 st
014-00129	Hudadapter (engångssensor 014-00321)	1 st
014-00322	Kärntemperatursensor, Vuxen, Engångs	20 st
014-00321	Hudtemperatursensor, Engångs	20 st
014-00020	Kärntemperatursensor, Vuxen, Flergångs	1 st
014-00021	Hudsensor, Flergångs	1 st
014-00005	Kärntemperatursensor, Barn, Flergångs	1 st
3500	Cure Wrap, Vuxen, en storlek	8 st
3518	Cure Wrap, Infant < 4 kg	8 st
3521	Cure Wrap, Infant 4-7 kg	8 st
SM008	Small 4 st. < 4 kg + 4 st. 4-7 kg	8 st
MD008	Medium 4 st. 7-11 kg + 4 st. 79-91 cm	8 st
LD008	Large 4 st. 91-104 cm + 4 st. 104-122 cm	8 st
XL008	X-large 4 st. 122-135 cm + 4 st. > 135 cm	8 st
3131	Universal ThermoWrap™ 79-91 cm	12 st
3136	Universal ThermoWrap™ 91-104 cm	12 st
3141	Universal ThermoWrap™ 104-122 cm	12 st
3148	Universal ThermoWrap™ 122-135 cm	12 st
3153	Universal ThermoWrap™ 135-152 cm	12 st
3160	Universal ThermoWrap™ 152-168 cm	12 st
3166	Universal ThermoWrap™ 168-180 cm	12 st
3363	ThermoWrap™, en storlek passar både operation och intensivvård	12 st